

# AMD Standard Diagnostic Video Scopes

All the tools you need...all connecting to a common platform

## AMD-2015 ENT Scope

A must for comprehensive ear, nose and throat (ENT) examinations, the AMD-2015 ENT/Otoscope combines the functionality of a high performance otoscope, a short sinus scope and an oral exam scope in a single diagnostic device to meet your ENT needs. Pneumatic Otoscopy for easy assessment of outer-ear canal, ear drum, and middle ear. Permits excellent imagery for lower sinus and upper throat examinations. Exceptionally clear, wide field of view. Connects to the AMD-400 and AMD-500 camera platforms. One of our most popular diagnostic devices, the ENT/otoscope is acclaimed for its use in otoscopy to cover the full range of ear, nose and throat examinations.



### Specifications

- **Optics:** Rod lens
- **Illumination:** POC Illumination lamp
- **Field of View:** Wide angle 80 degrees
- **Dimensions:** 2.7mm diameter, 80mm in length
- **Specula:** Adult uses 4mm specula; pediatric uses 2.5mm specula

## AMD-2030 Dermoscope

A superior skin surface contact microscope for use in primary care, pediatrics or dermatology. Easy, integrated focus controls. Removable contact plate provides easy cleaning and sterilization and eliminates cross-contamination. Millimeter scale in contact plate allows lesion measurement. Connects to the AMD-400 and AMD-500 camera platforms. The AMD-2030 is widely used in primary care and pediatrics, and by dermatology specialists, to diagnose dermatologist conditions and to document and archive lesions.



### Specifications

- **Optics:** Integrated focusing optics with adjustable focus
- **Illumination:** POC Illumination lamp
- **Field of View:** 13mm