# AMD-2040 Nasopharyngoscope



## Features

 An internal insertion tube strain relief that increases the working length of the Nasopharyngoscope, while still providing protection of the upper section of the insertion tube

•The AMD-2040 connects to the AMD-400s and AMD-500 Camera & Illumination System. Also requires the AMD-2451 23mm Coupler.

# **Specifications**

- 3.8mm Outer Diameter:
- •Depth of Field: 5mm-Infinite
- •Working Length:
- •Tip deflection: up/down
- •Field of View:

25mm

125°

Wide angle

The AMD-2040 Nasopharyngoscope is designed for examination of the upper airway from the nasal passage to the vocal cords.

#### **Insertion Tube:**

•Depth Markers: Depth marks are located at 10cm and 20cm from the distal tip of the insertion tube •Bending Section: The bending


section deflects up and down when the articulation lever is rotated

•Distal Tip: The termination of the image fiber bundle (at the objective lens) & the illumination fiber bundles (light guides)

## **Control Body:**

•Eyepiece: Allows the user to connect to the AMD-400s and AMD-500 Camera and Illumination System

•Focus Ring: The focus ring adjusts the image to the user's eyesight







Rev 1.1